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CE

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GoMed
PARTNER  medizintechnik

GE KISS Leadwires Instructions For Use

Functional Description General Information

The GE Healthcare Electrode Application System KISS™ is used to assist in the application of ECG electrodes for the acquisition of resting and exercise ECGs. The system is optimally suited from an ergonomic viewpoint.

The electrode leads are conducted to a common trunk cable via a central electrode distributor. Integrated in the electrodes - and trunk cable - is the pneumatic tubing connector.

A number of instrument versions are available to match the ECG recorder and acquisition module used.

An option for the electrode application system KISS™ comprises a pole with an extendible swivel arm adjustable in length and height for installation on a GE Healthcare instrument cart. With the brackets provided in the accessory kit, the system can be mounted to a wall or on an instrument table.

An electrically operated suction pump, installed in the ECG recorder / stress system or in a system power supply unit generates a constant negative pressure, the negative pressure controller in the distributor enabling controlled adjustment.


The negative pressure can be adjusted within a range of 80 to 220 mbar in accordance with the required application.

The suction pump only operates when negative pressure must be created after attaching the electrodes or when a loss in pressure for electrodes already attached has to be compensated for. This permits the use of a quiet, low Watt rating pump.

The power supply unit supplies the suction pump with mains voltages from 110V to 240 V.

The silver/silver chloride electrodes are attached by suction by pressing briefly on the suction cup at the application site. A valve opens in the electrode lead, enabling the generation of negative pressure.

Operating Controls and System Components

Application System KISS and pole

Item	Description	Item	Description
1	Electrode distributor	11	Power connector for pump module
2	Negative pressure control	12	Pump module
3	Electrode leads with valve	13	Screw clamp for fixation of pole to table top
4	Silver silver chloride electrode	14	Wall-mount facility for pole
5	Filter disk => electrode	15	Pole top with retaining screw, accommodates contact spray bottle
6	Suction cup => electrode	16	Locking catch to adjust the height of the swivel arm
7	Electrocardiograph and pump connecting cable (trunk cable)	17	Locking catch to adjust the length of the swivel arm
8	Connector for electrocardiograph	18	Bracket for electrode distributor
9	Luer-Lock connector for suction pump	19	Electrode lead hanger
10	Extension tube	20	Power cord hanger

Cleaning and Disinfecting the Patient and Leadwires

Disconnect the cables from the system before cleaning and disinfecting them.

Clean the cables with a cloth moistened with soap water. Use a disinfectant for disinfection of the leadwires and cables.

CAUTION:
EQUIPMENT DAMAGE - Under no circumstances should leadwires be immersed in the cleaning solution or disinfectant. Even if the system were not under negative pressure, the liquid would enter the system, spoiling the filter disks and clogging the pneumatic system.

Cleaning Electrodes and Electrode Cable

For a thorough cleaning the electrode cable (3) can be disconnected from the electrode (4). When doing this it is practical to disconnect the cable plug.

The suction cup (6) can now be removed from the electrode. This permits access to the filter disk (5), which can be replaced.

The surface of electrodes with minor contamination can be cleaned with a cloth moistened with water and detergent.

CAUTION:
EQUIPMENT DAMAGE - Do not scrape electrodes with metal objects. The silver-silver chloride structure could be damaged by this and lose its stabilizing effect on the electrode interface.

To disinfect the electrodes, wipe them down with a cloth moistened with 80% alcohol. All disinfectants containing alcohol can be used for disinfection.

CAUTION:

EQUIPMENT DAMAGE - After use do not immerse the electrodes into the cleaning solution. The liquid entering the negative pressure system would soak and hence spoil the filter disk.

For a system in frequent use replace the filter disks once a week!

Replacing Leadwires

The leadwires are subject to wear and tear and, therefore, they can be replaced. The leadwires are available without the electrodes as spare parts.

The leadwires are plugged into the distributor head for easy replacement. When removing the lead pull on the plug. When plugging in a new lead do not use force; the pneumatic and electrical plug-in connectors could become damaged.

Spare Parts List

Single Leadwires for KISS and KISS Multilead

X-GEK-127SC-C1	Leadwire C1
X-GEK-127SC-C2	Leadwire C2
X-GEK-127SC-C3	Leadwire C3
X-GEK-127SC-C4	Leadwire C4
X-GEK-127SC-C5	Leadwire C5
X-GEK-127SC-C6	Leadwire C6
X-GEK-147SC-F	Leadwire F
X-GEK-147SC-L	Leadwire L
X-GEK-147SC-R	Leadwire R
X-GEK-147SC-N	Leadwire N
GEK10-AD-I	Trunk cable for 10-lead KISS, standard
Y05-T286	Suction electrode for GE KISS
	Filter disk

Environment

Operation

- temperature between +10 to +40 °C
- relative humidity between 25 and 95%, non-condensing
- atmospheric pressure between 700 and 1060 hPa

Storage and transport

- KISS w/o pump:
Temperature between +5°C and +45°C
- relative humidity between 10 and 95%, non-condensing
- atmospheric pressure between 500 and 1060 hPa

Weight

- electrode distributor with 10 electrodes approx. 800 g



Negative pressure adjustment for pump module continuously adjustable between 80 and 220 mbar; in position

⊖ negative pressure < 30 mbar.

length of electrode leads limbs 140 cm chest 120 cm.

length of negative pressure tube connected to separate pump or ECG Device: 270 cm.

Symbol explanation

Caution	Lot number	Manufacturer	Latex free	See the instructions
				
Waste Electrical and Electronic Equipment	Date of manufacture	Catalogue number	Authorized representative in the European Community	Keep away from sunlight
				
Keep dry	Minimum/Maximum Temperature			
				

Unimed Medical Supplies, Inc

Bld#8, Nangang 3rd Industrial Park, Tangtuo, Shiyuan, 518108
ShenZhen, People's Republic of China

Tel: (86-755) 26695165

Fax: +86 755 26697984

Web: www.unimed.cn

Email: info@unimed.cn

REF: Y07-M197
Rev.: A0

Obelis s.a.
Boulevard Général Wahis 53,
1030 Brussels, BELGIUM
Tel: +(32) 2. 732.59.54
Fax: +(32) 2.732.60.03
E-Mail: mail@obelis.net

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